



**Reading Jug**


Write the reading shown by each jug.




\_\_\_\_\_ mL




\_\_\_\_\_ mL




\_\_\_\_\_ mL




\_\_\_\_\_ mL




\_\_\_\_\_ mL




\_\_\_\_\_ mL



\_\_\_\_\_ mL












\_\_\_\_\_ mL



\_\_\_\_\_ mL

**Answer key**

Write the reading shown by each jug.

350 mL700 mL250 mL650 mL450 mL900 mL800 mL100 mL550 mL